

Gal-G2-CNP



成为世界一流的特种酶产品与服务提供商
info@hzymes.com
www.hzymes.com

Manual Version: V1.0

Product name

Gal-G2-CNP
(2-Chloro-4-nitrophenyl-4-O- β -
-D-galactopyranosylmaltoside)

Abbreviation

Gal-G2-CNP

Cat. No.

HH8001


CAS

157381-11-8

Appearance


Off-White to pale yellow powder

Constitutional formula



Product function

The substance is α - Starch hydrolase substrate



Molecular formula / molecular weight

C₂₄H₃₄ClNO₁₈ / 659.97

Purity

≥90 %

Detection wavelength

λ max (Color wavelength) = 405 nm

Shipping and storage

-25 ~ -15 °C, sealed, dry, light avoidance storage.

Storage stability

2 years

Product advantages

- 1) Good water solubility;
- 2) High sensitivity;
- 3) Good linearity;
- 4) Good stability

Precaution

Only for in vitro diagnostic products and shall not be used in human or in vivo experiments.